I-Series Peritoneal Catheter  
Patient Information Packet

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### I-Series Peritoneal Catheter Patient Information Packet

#### PIP Statement:

This patient information packet is intended to provide general information regarding this device and should not replace the Information for Use.

#### Information for Users/Healthcare Professionals:

The following information is intended for users/healthcare professionals. Following that information is a summary intended for patients.

#### Indications for Use:

The Jet I-Series Peritoneal Dialysis Catheter is indicated for short-term or long-term use where peritoneal access is required for the purpose of peritoneal dialysis for Acute Kidney Injury (AKI) or Chronic Kidney Disease (CKD).

#### Device Lifetime:

Medcomp® catheters are subjected to, and must pass, simulated use testing intended to replicate 12 months use as part of device development. The I-Series Peritoneal Dialysis Catheter passed this testing. Although Medcomp® catheters materials contain non-degradable polymers, fully functional catheters may be removed for other reasons,

such as intractable infection or change of therapy. Published clinical literature does not always focus on the physical lifetime of a catheter for these reasons. In the case of the Jet I-Series Peritoneal Dialysis Catheter, 131 catheters had a 20.71 month [95% CI: 17.76 - 23.66 months] duration of use that has been found in clinical use reported to date.<sup>1</sup>

Based on this information, the Jet I-Series Peritoneal Dialysis Catheter has a 12-month lifetime; however, the decision to remove and/or replace the catheter should be based on clinical performance and need, and not any predetermined point in time.

#### MR Safety Information:

 MR Conditional – 1.5 Tesla and 3 Tesla (artifacts may present imaging problems if MRI area of interest is on or near area where device is located).

Report Conclusion: MRI Information MR Conditional.

Non-clinical testing has demonstrated that the peritoneal catheter system is MR Conditional. A patient with this device can be safely scanned in an MR system with the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0 Tesla (3.0T).
- Maximum spatial gradient field of 19 T/m (1900 G/cm).
- Maximum MR system reported, whole-body averaged specific absorption rate (SAR) of 2.0 J/kg (normal operating mode).

Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than or equal to 2.7 °C after 15 minutes of continuous scanning. In non-clinical testing, the image artifact caused by the device extends up to 6.1 cm when imaged with a gradient echo pulse sequence in a 3.0T MR system.

#### Site Care:

Clean skin around catheter. Chlorhexidine gluconate solutions are recommended; however, iodine-based solutions can also be used. Cover the exit site with occlusive dressing and leave extensions, clamps, and caps exposed for access by staff. Wound dressings must be kept clean and dry. **Patients must not swim, shower, or soak dressing while bathing.** If profuse perspiration or accidental wetting compromises adhesion of dressing, the medical or nursing staff must change the dressing under sterile conditions.

- Catheter is compatible with ointments.
- Clean skin around catheter. Chlorhexidine gluconate solutions are recommended.

Alternate compatible solutions/ointments include:

- Chloraprep<sup>®</sup>, Chlorhexidine Gluconate 2%
- Isopropanol 70%
- Bactroban<sup>®</sup>, Mupirocin Ointment 2%
- Silvadene<sup>®</sup>, 1% Silver Sulfadiazine Cream

Solutions should be allowed to completely dry before applying an occlusive dressing.

### Warnings, Precautions, or Measures to be Taken by Healthcare Professional:

- Before attempting procedure, ensure that you are familiar with the potential complications and their emergency treatment should any of them occur.
- Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care professionals should always use Universal Blood and Body Fluid Precautions in the care of all patients.
- The medical techniques and procedures described in these instructions do not represent all medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.
- Call 215-256-4201 or visit [www.medcomp.net](http://www.medcomp.net) for any information to ensure safe use of the device, warnings, precautions or measures to be taken by the patient or a healthcare professional.

### How To Fill Out Patient ID Card: Instructions for Completion of Patient ID Card:

1. Name of the patient or patient ID. To be filled by the healthcare institution/provider.
2. Date of implantation. To be filled by the healthcare institution/provider.
3. Name and address of the healthcare institution/provider. To be filled by the healthcare institution/provider.

 [www.medcomp.net/patientinformation](http://www.medcomp.net/patientinformation)

 \_\_\_\_\_

 \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Medical Components, Inc.  
14901 Duple Dr., Independence, PA 15118, U.S.A.  
P: 215 256 4201 F: 215 256 1787  
[www.medcomp.net](http://www.medcomp.net)

### PATIENT INFORMATION:

The information presented below is intended for patients or lay persons. A more extensive summary of this device prepared for healthcare professionals is found in the first part of this document.

Contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

### Site Care:

Clean skin around catheter. Chlorhexidine gluconate solutions are recommended; however, iodine-based solutions can also be used. Cover the exit site with occlusive dressing and leave extensions, clamps, and caps exposed for access by staff. Wound dressings must be kept clean and dry. **Patients must not swim, shower, or soak dressing while bathing.** If profuse perspiration or accidental wetting compromises adhesion of dressing, the medical or nursing staff must change the dressing under sterile conditions.

- Catheter is compatible with ointments.
- Clean skin around catheter. Chlorhexidine gluconate solutions are recommended.

Alternate compatible solutions/ointments include:

- Chloraprep<sup>®</sup>, Chlorhexidine Gluconate 2%
- Isopropanol 70%
- Bactroban<sup>®</sup>, Mupirocin Ointment 2%
- Silvadene<sup>®</sup>, 1% Silver Sulfadiazine Cream

Solutions should be allowed to completely dry before applying an occlusive dressing.

### Warnings, Precautions, or Measures to be Taken by Patient:

- To reduce the risk of bacterial contamination of catheter, wear a mask over your nose and mouth and wash your hands before touching your catheter. Observe proper sterile technique.
- Never put anything on your exit site without direction from a healthcare professional.
- Ask the doctor to explain the signs and symptoms of infection.
- Call a healthcare professional immediately if you think you have signs of an infection including: redness, swelling, pain or tenderness, drainage.

[www.medcomp.net/patientinformation](http://www.medcomp.net/patientinformation)

### Materials or Substances the Patient May be Exposed to:

The percentage ranges in the table below are based on the weights of the 18 cm catheter (4.34 g) and the 63 cm catheter (12.40 g).









Material	% Weight (w/w)
Silicone	56.65 - 78.20
Barium Sulfate	11.28 - 13.80
Nylon	7.33 - 20.94
Polyester	0.35 - 11.13

**Note:** Stainless steel containing accessories such as tunneler, stylet, introducer needle, guidewire and scalpel may contain up to 1.0% weight of the CMR substance cobalt. This does not apply to Catheter Only sets.

### References:

1. Medcomp. PMCF\_PD\_202\_Report\_B. 20 December 2021.

### SYMBOL TABLE

	Manufacturer *
	Serial Number *
	Batch Code *
	Authorized Representative in the European Community/European Union*
	Medical Device *
	Unique Device Identifier *
	MR Conditional - 1.5, 3.0 Tesla ****
	Patient Information Website *
	Date *
	Healthcare Center or Doctor *
	Patient Identification *
	United Kingdom Responsible Person **
	Consult Instructions for Use or Consult Electronic Instructions for Use *
	Catalogue Number *
	Date of Manufacture *
	Customer Feedback **
	Contains Hazardous Substances *
	Contains the CMR Substance, Cobalt

\*This symbol is in accordance with ISO 15223-1.

\*\*Not a recognized symbol.

\*\*\*This symbol is in accordance with ASTM F 2503.

Medcomp<sup>®</sup> is a trademark of Medical Components, Inc. registered in the United States.

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The diagram shows a patient information card with three main sections. The first section is for the patient's name, with a question mark icon and a '31' icon. The second section is for the date of insertion, with a calendar icon. The third section is for the hospital name, with a person icon. The URL [www.medcomp.net/patientinformation](http://www.medcomp.net/patientinformation) is printed at the top right of the card.

Medical Components, Inc.  
1699 Dale Dr., Hazelville, IA 50438, U.S.A.  
Tel.: 215-256-4000 • Fax: 215-256-1787  
[www.medcomp.net](http://www.medcomp.net)

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© 2015 256-4201 P | 215-256-1787
Medcomp®, Medical Components, Inc.'in ABD'deki tescilli ticari markasıdır.
Harleysville, PA 19438, U.S.A.
MPS Medical-Product-Service UK Ltd
MPS Medical Product Service GmbH