

Jet Medical SA now offers straight and curled Peritoneal Dialysis Catheters in both I-Series and V-Series (pre-curved) configurations. These catheters are manufactured using only the highest quality medical-grade silicone and are available with your choice of a single subcutaneous cuff or double intraperitoneal cuffs. The large, strategically-placed side holes around each catheter's intraperitoneal section facilitate increased dialysis solution flows and a reduction in the incidence of omental occlusion. Our special Jet Peel Cath Sheath Introducer is supplied as a separate set, and allows these catheters to be placed percutaneously- without the need for an elaborate surgical procedure or the use of trocars. Accessories include: caps, stylets, clamps, and adaptors, as well as the 16FR Peel Cath Introducer Set.



Features & Benefits

- Available in straight or curled configurations.
- Available with single or double pre-attached dacron cuffs.
- Manufactured from medical-grade silicone.
- A highly radiopaque stripe runs the full length of each catheter, improving radiographs and catheter placement confirmation.
- Large side holes placed strategically around each catheter's intraperitoneal section provide excellent dialysis solution flows and lower incidence of omental occlusion.
- Long external lumen makes catheter maintenance simple.
- Peel Cath Sheath Introducer Set eliminates the need for surgical implantation or the use of trocars.

ORDERING INFO::

Catheters Contain:

- (1) Catheter
- (1) Clamp
- (1) Female Luer
- (1) Male Luer End Cap

Sets Contain:

- (1) Catheter
- (1) 18ga Introducer Needle
- (1) 10cc Syringe
- (1) #11 Scalpel
- (1) 16F Sheath/Dilator
- (1) Tunneling Tool
- (6) Gauze Sponges
- (1) .038 Guidewire
- (1) Beta Clamp
- (1) Adaptor
- (1) Male Luer

CATHETERS 5/BOX

SETS 5/BOX

Jet Peritoneal Catheters

I-SERIES

JDC31S		31cm Straight Catheter with Double Cuff
JDC32S		32.25cm Straight Catheter with Double Cuff
JDC37S		37cm Straight Catheter with Double Cuff
JDC42S	JET-242	42cm Straight Catheter with Double Cuff
JDC47S		47cm Straight Catheter with Double Cuff

JDC57C	JET-257	57.5cm Curled Catheter with Double Cuff
JDC63C	JET-263	63cm Curled Catheter with Double Cuff

JSC31S		31cm Straight Catheter with Subcutaneous Single Cuff
JSC37S		37cm Straight Catheter with Subcutaneous Single Cuff
JSC42S		42cm Straight Catheter with Subcutaneous Single Cuff
JSC46S	JET-146S	46cm Straight Catheter with Subcutaneous Single Cuff

JPP41S	JET-141	41cm Straight Catheter with Preperitoneal Single Cuff
JPP46S	JET-146	46cm Straight Catheter with Preperitoneal Single Cuff

JPP40C		40.25cm Coiled Catheter with Preperitoneal Single Cuff
JPP57C		57.5cm Coiled Catheter with Preperitoneal Single Cuff
JPP60C	JET-160	60cm Coiled Catheter with Preperitoneal Single Cuff

V-SERIES

JS200101	JS200301	Straight, 43cm Left V-NECK PD Catheter with Double Cuff
JS200102	JS200302	Straight, 43cm Right V-NECK PD Catheter with Double Cuff
JS200201	JS200401	Coiled, 63cm Left V-NECK PD Catheter with Double Cuff
JS200202	JS200402	Coiled, 63cm Right V-NECK PD Catheter with Double Cuff

Jet Peel Cath Introducer Set Contains:

- (1) Peel Cath Introducer
- (1) Guidewire
- (1) Introducer Needle
- (2) 10cc Syringes

PART NO.

ACCESSORIES

JPC216	16F Jet Peel Cath Introducer Set
JMEC	Male Luer Lock End Cap
JBCA	Barbed Luer Lock Catheter Adaptor
JBC	Beta Clamp
JCS	Catheter Stylet
JTS	Tunneling Stylet
JTSR	Tunneling Stylet (Rigid)
JTC	Titanium Connector (male-female)

Jet
Medical^{SA}

Rue Fritz-Courvoisier 40
2300 La Chaux-de-Fonds
Switzerland
Tel: +41-32-9679021
Fax: +41-32-9679001
www.jetmedical.net
customerservice@jetmedical.net