



JET INFUSION

The Jet Infusion catheter allows the administration of nutritional replacement therapy and the infusion of medications and solutions. Double, triple, and quad lumens allow multiple drugs to be administered through the same catheter.

FEATURES & BENEFITS

Adjustable Wing allows securement at puncture site.

Printed depth markings.

Vessel dilators aid in placing the catheter percutaneously.

CVC's are available in a variety of lengths, French sizes, and diameters.

ORDERING INFORMATION

JET INFUSION CATHETERS (2,3,4 LUMEN)

Catalog #	Description	Qty
JETIC245	4FR x 5cm Double Lumen Catheter	10/Box
JETIC248	4FR x 8cm Double Lumen Catheter	10/Box
JETIC2413	4FR x 13cm Double Lumen Catheter	10/Box
JETIC2430	4FR x 30cm Double Lumen Catheter	10/Box
JETIC258	5FR x 8cm Double Lumen Catheter	10/Box
JETIC2513	5FR x 13cm Double Lumen Catheter	10/Box
JETIC2530	5FR x 30cm Double Lumen Catheter	10/Box
JETIC358	5.5FR x 8cm Triple Lumen Catheter	10/Box
JETIC355	5.5FR x 13cm Triple Lumen Catheter	10/Box
JETIC3530	5.5FR x 30cm Triple Lumen Catheter	10/Box
JETIC2715	7FR x 15cm Double Lumen Catheter	10/Box
JETIC2720	7FR x 20cm Double Lumen Catheter	10/Box
JETIC2730	7FR x 30cm Double Lumen Catheter	10/Box
JETIC3715	7FR x 15cm Triple Lumen Catheter	10/Box
JETIC3720	7FR x 20cm Triple Lumen Catheter	10/Box
JETIC3730	7FR x 30cm Triple Lumen Catheter	10/Box

MULTI-LUMEN SETS CONTAIN:

- (1) Catheter with Clamps
- (1) J/Flex Guidewire
- (1) Introducer Needle
- (2) End Caps
(3 for triple lumen, 4 for quad lumen)
- (1) 5cc Syringe
- (1) Vessel Dilator
- (1) Movable Wing

SINGLE-LUMEN SETS CONTAIN:

- (1) Catheter with Clamps
- (1) J/Flex Guidewire
- (1) Introducer Needle
- (1) Vessel Dilator
- (1) Movable Wing

JET INFUSION CATHETERS (1 LUMEN)

Catalog #	Description	Qty
JET11415	14GA x 15cm Single Lumen Catheter	20/Box
JET11420	14GA x 20cm Single Lumen Catheter	20/Box
JET11430	14GA x 30cm Single Lumen Catheter	20/Box
JET11615	16GA x 15cm Single Lumen Catheter	20/Box
JET11620	16GA x 20cm Single Lumen Catheter	20/Box
JET11630	16GA x 30cm Single Lumen Catheter	20/Box
JET11815	18GA x 15cm Single Lumen Catheter	20/Box
JET11820	18GA x 20cm Single Lumen Catheter	20/Box
JET11830	18GA x 30cm Single Lumen Catheter	20/Box
JET12013	20GA x 13cm Single Lumen Catheter	20/Box
JET12413	24GA x 14cm Single Lumen Catheter	20/Box