



## JET INFUSION

**Central Vein Catheters facilitate infusion therapy** in the critical care environment by permitting the infusion of various drugs and antibiotics, blood sampling, and central venous monitoring. Double and triple lumen models allow multiple drugs to be administered through the same catheter.

The Jet Infusion Central Vein Catheters provide a vital tool for the intensivist where central access is limited but vital. A small guidewire (.018") and needle (21g) facilitate safe, small vein entry with the smaller bore catheters.

Exclusive techniques in resin fabrication have brought together high flow lumens and soft tips into several multi-lumen configurations (double 4FR, 5FR and triple 5.5FR). Soft polyurethane lumens with super soft tips are also available in smaller French size configurations:

- 4FR Double Lumen with 20/19g Lumens
- 5FR Double Lumen with 19/18g Lumens
- 5.5FR Triple Lumen with 22/22/20g Lumens

## FEATURES & BENEFITS

Movable wing/clamp allows anchorage at puncture site regardless of depth of catheter. This minimizes trauma and irritation to the exit site.

Printed depth markings assist in accurate central venous placement from the right or left subclavian or jugular vein approaches.

Atraumatic tips are less traumatic to the vessel, minimizing vessel erosion, hemothorax, and cardiac tamponade. Available in double, triple, and quad lumen versions.

Radiopacity allows fluoroscopic confirmation of catheter placement. The multi-lumen versions are more radiopaque at the tip to easily confirm fluoroscopic tip placement.

Vessel dilators allow "super soft" catheters to be easily placed percutaneously, thereby ensuring a soft catheter/vein interface.

Quality & Economy - sold at competitive prices and manufactured under world class requirements.

- FDA/GMP Compliant
- ISO13485 Certified Quality System
- CE (European) Certified
- Japanese Shonin

# ORDERING INFORMATION

## JET INFUSION CATHETERS (2,3,4 LUMEN)

CATALOG NO. (10/BOX)	DESCRIPTION
JETIC245	4FR x 5cm Double Lumen Catheter
JETIC248	4FR x 8cm Double Lumen Catheter
JETIC2413	4FR x 13cm Double Lumen Catheter
JETIC2430	4FR x 30cm Double Lumen Catheter
JETIC258	5FR x 8cm Double Lumen Catheter
JETIC2513	5FR x 13cm Double Lumen Catheter
JETIC2530	5FR x 30cm Double Lumen Catheter
JETIC358	5.5FR x 8cm Triple Lumen Catheter
JETIC355	5.5FR x 13cm Triple Lumen Catheter
JETIC3530	5.5FR x 30cm Triple Lumen Catheter
JETIC2715	7FR x 15cm Double Lumen Catheter
JETIC2720	7FR x 20cm Double Lumen Catheter
JETIC2730	7FR x 30cm Double Lumen Catheter
JETIC3715	7FR x 15cm Triple Lumen Catheter
JETIC3720	7FR x 20cm Triple Lumen Catheter
JETIC3730	7FR x 30cm Triple Lumen Catheter
JETIC2816	8FR x 16cm Double Lumen Catheter
JETIC2820	8FR x 20cm Double Lumen Catheter
JETIC4820	8.5FR x 20cm Quad Lumen Catheter

### MULTI-LUMEN SETS CONTAIN:

- (1) Catheter with Clamps
- (1) J/Flex Guidewire
- (1) Introducer Needle
- (2) Injection Caps  
(3 for triple lumen, 4 for quad lumen)
- (1) 5cc Syringe
- (1) Vessel Dilator
- (1) Movable Wing
- (1) Wing Clamp

### SINGLE-LUMEN SETS CONTAIN:

- (1) Catheter with Clamps
- (1) J/Flex Guidewire
- (1) Introducer Needle
- (1) Vessel Dilator
- (1) Movable Wing
- (1) Wing Clamp

**Jet**  
Medical<sup>®</sup> SA

Rue Fritz-Courvoisier40  
2300 La Chaux-de-Fonds  
Switzerland  
Tel: +41-32-9679021  
Fax: +41-32-\*9679001  
www.jetmedical.net  
customerservice@jetmedical.net

## JET INFUSION CATHETERS (1 LUMEN)

CATALOG NO. (20/BOX)	DESCRIPTION
JET11415	14GA x 15cm Single Lumen Catheter
JET11420	14GA x 20cm Single Lumen Catheter
JET11430	14GA x 30cm Single Lumen Catheter
JET11615	16GA x 15cm Single Lumen Catheter
JET11620	16GA x 20cm Single Lumen Catheter
JET11630	16GA x 30cm Single Lumen Catheter
JET11815	18GA x 15cm Single Lumen Catheter
JET11820	18GA x 20cm Single Lumen Catheter
JET11830	18GA x 30cm Single Lumen Catheter
JET12013	20GA x 13cm Single Lumen Catheter
JET12413	24GA x 13cm Single Lumen Catheter