




LONG-TERM HEMODIALYSIS PRODUCT CATALOG



Jet Medical SA develops, manufactures, markets, and supports cutting-edge vascular access devices and accessories to meet the clinical needs of the medical industry. Our Long-Term Hemodialysis product lines offer a variety of options for patients. This includes single and dual lumen catheters for antegrade or retrograde insertion techniques. These catheters are available in a wide range of materials, tip designs, French sizes, lengths, and kit configurations.

Message to the Customers

Jet Medical SA's vision is to be the leading solution provider for every customer we serve. Our mission, to drive sustainable growth of our business through the continued investment in our team members, technology and quality while delivering world-class products and services. Since the founding of Jet Medical SA we have built pioneering medical devices in the field of Dialysis. Continued engagement between our in-house clinical staff and partnerships with top physicians in the field of nephrology and vascular access drive innovation of new device designs to serve our customers every need.



www.JetMedical.net

LONG TERM HEMODIALYSIS

DUO-JET® II
JET FLOW
REPAIR KIT

IMPORTANT RISK INFORMATION

DUO-JET® II

Indications for Use: The Duo-Jet® II Catheter is indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis.

Contraindications: Known or suspected allergies to any of the components of the catheter or the kit. This device is contraindicated for patients exhibiting severe, uncontrolled coagulopathy or thrombocytopenia.

JET FLOW

Indications for Use: The Jet Flow Catheter is indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis.

Contraindications: Known or suspected allergies to any of the components of the catheter or the kit. This device is contraindicated for patients exhibiting severe, uncontrolled coagulopathy or thrombocytopenia.

JET MEDICAL REPAIR KIT

Indications for Use: The Jet Medical Long-Term Catheter Repair Kit is indicated to replace catheter luer connectors, clamps, or repair extensions.

Contraindications: Known or suspected allergies to any of the components of the repair kit.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation




SINGLE ROUND INTERNAL
LUMEN DESIGN

Polyurethane	Single Round	10F	Straight	Retrograde
Material	Tip Design	French Size	Configuration	Insertion Type

DJLT2000L	10F x 52CM DUO-JET [®] II CATHETER, LEFT SIDE APPROACH	5/BOX
DJLT2000R	10F x 52CM DUO-JET [®] II CATHETER, RIGHT SIDE APPROACH	5/BOX

Set Contents: (2) Catheter (2) Introducer Needle (2) Guidewire w/ Advancer (2) Tunneler
(1) Coaxial Introducer Sheath (2) Peelable Introducer (1) Arterial Extension Set (1) Venous
Extension Set (2) Clamp (2) Catheter Plug (2) End Cap (1) Catheter Securement Device
(1) Patient ID Card (1) Patient Information Packet



Refer to the Table of Contents for Important Risk Information regarding this device.

Polyurethane	Step Tip	14.5F	Straight, Pre-Curved	Antegrade
Material	Tip Design	French Size(s)	Configurations	Insertion Type




DOUBLE "D" INTERNAL
LUMEN DESIGN

14.5F CATHETER SET

JFC1424	14.5F x 24CM JET-FLOW CATHETER, STRAIGHT	5/BOX
JFC1428	14.5F x 28CM JET-FLOW CATHETER, STRAIGHT	5/BOX
JFC1432	14.5F x 32CM JET-FLOW CATHETER, STRAIGHT	5/BOX
JFC1436	14.5F x 36CM JET-FLOW CATHETER, STRAIGHT	5/BOX
JFC1440	14.5F x 40CM JET-FLOW CATHETER, STRAIGHT	5/BOX
JFC1424PC	14.5F x 24CM JET-FLOW CATHETER, PRE-CURVED	5/BOX
JFC1428PC	14.5F x 28CM JET-FLOW CATHETER, PRE-CURVED	5/BOX
JFC1432PC	14.5F x 32CM JET-FLOW CATHETER, PRE-CURVED	5/BOX

Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire w/ Advancer
(1) Tunneler (1) Dilator (1) Valved Peelable Introducer (2) End Cap (1) Patient ID Card
(1) Patient Information Packet



Refer to the Table of Contents for Important Risk Information regarding this device.

ARA101	JET MEDICAL REPAIR KIT	5/BOX
--------	------------------------	-------

Kit Contents: (1) Luer Assembly (1) Scissors (1) Adhesive Backed Drape (3) Clamp (2) End Cap

*Repairs JET-FLOW

ORDERING INFO & POLICIES

All orders should be directed to Jet Medical's Customer service department. Orders should be for full box quantities, accompanied with purchase order number.

Products are shipped via United Parcel Service or Fed Ex and are F.O.B. Origin in the Continental United States. Express service is also available at the expense of the customer.

Export shipments are shipped air freight to maintain traceability. Should a customer wish air parcel post shipments, Jet Medical SA will not assume responsibility for lost shipments. All shipments outside of the Continental United States are EX Works Origin.

Terms

All orders will be billed at prevailing prices at the time of shipment.

Terms are net 30 days for domestic accounts. Outside the U.S., terms are prepayment by wired funds or letter of credit at sight.

Prices are in U.S. dollars, and are subject to tax, if any. All prices are subject to change without notice.

Return Goods Policy

Return goods must be authorized by Jet Medical SA prior to return. Jet Medical SA will accept return of any standard product within 30 days, upon approval, providing it is undamaged and in original sealed package, and returned prepaid. Custom products are nonreturnable.

Warranty

- Medcomp® warrants its products to be free from defects in workmanship and materials for the period up to and including the expiration date.
- Medcomp® shall not be liable for any incidental, special or consequential loss, damage or expense directly or indirectly arising from the use of its products.
- Liability is limited to replacing, at Medcomp's option, any products which shall under normal use and service appear to the company to have been defective in materials or workmanship.
- Damage to the products through misuse, neglect, or accident voids this warranty.
- No agent, employee or representative of Medcomp® has the authority to bind Medcomp® to any other warranty.
- This warranty is expressly in lieu of any other implied warranties, including any merchantability of fitness and of any other obligation on the part of the seller.
- NOTE: Because of continuing product improvement, prices, specifications, and model availability are subject to change without notice.
- Medcomp® is a trademark of Medical Components, Inc. registered in the United States.
- Duo-Jet® is a trademark of Jet Medical S.A. registered in the United States.

Patents

- To view a full list of our patents, please visit www.medcomp.net/patent/



LONG-TERM HEMODIALYSIS PRODUCT CATALOG



Rue Fritz-Courvoisier 40
2300 La Chaux-de-Fonds
Switzerland
Tel: +41-32-9679021
Fax: +41-32-9679001
www.jetmedical.net
customerservice@jetmedical.net



Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201 F: 215.256.1787
www.medcomp.net

ISO13485
Certified Quality System

PN 2685LTHD REV. A 10/22