

JET LT CVC

Central Venous Catheter



FEATURES & BENEFITS

Lumens

Single, Dual & Triple

Lengths

90cm

Material

Silicone

Available In

7F, 9.5F, 11F, & 12.5F


JET LT CVC

Central Venous Catheter

| JET LT CVC BASIC SET | | |
|----------------------|---------------------------|--------|
| JETHC17 | 7F X 90CM SINGLE LUMEN | 10/BOX |
| JETHC19 | 9.5F X 90CM SINGLE LUMEN | 10/BOX |
| JETHC92 | 9.5F X 90CM DUAL LUMEN | 10/BOX |
| JETHC112 | 11F X 90CM DUAL LUMEN | 10/BOX |
| JETHC123 | 12.5F X 90CM TRIPLE LUMEN | 10/BOX |

Basic Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire
(1) Peelable Introducer* (1) Tunneler (1) Scalpel (1|2|3) End Cap(s)

*7F kit has a 8F Peelable Introducer
9.5F kit has a 10F Peelable Introducer
11F kit has a 12F Peelable Introducer
12.5F kit has a 13F Peelable Introducer



IMPORTANT RISK INFORMATION

Indications For Use: The Jet LT CVC Central Venous Catheter is designed for long-term central venous access. It can be used for total parenteral nutrition (TPN), infusion of I.V. liquids, blood, blood products and drugs. It can also be used for repeated withdrawal of blood samples.

Contraindications: This catheter is **NOT** intended for any use other than that which is indicated. It is **NOT** recommended for use in hemodialysis or in hemoperfusion procedures, nor in patients with severe chronic obstructive lung disease.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications. For patient information please visit www.medcompnet.com/patientinformation