

# JET FREE FLOW® ST

SHORT-TERM HEMODIALYSIS CATHETER



## FEATURES & BENEFITS

**Internal Lumen Design**  
Coaxial

**Material**  
Polyurethane

**Tip Design**  
Coaxial

**Available In**  
12.5F, 15.5F

**Configurations**  
Straight & IJ

# JET FREE FLOW® ST

SHORT-TERM HEMODIALYSIS CATHETER

12.5F Basic Set		
JFFS1215S	12.5F X 15CM JET FREE FLOW® ST, STRAIGHT	5/BOX
JFFS1220S	12.5F X 20CM JET FREE FLOW® ST, STRAIGHT	5/BOX
JFFS1224S	12.5F X 24CM JET FREE FLOW® ST, STRAIGHT	5/BOX
JFFS1212IJ	12.5F X 12CM JET FREE FLOW® ST, IJ	5/BOX
JFFS1215IJ	12.5F X 15CM JET FREE FLOW® ST, IJ	5/BOX
JFFS1220IJ	12.5F X 20CM JET FREE FLOW® ST, IJ	5/BOX

15.5F Basic Set		
JFFS1515S	15.5F X 15CM JET FREE FLOW® ST, STRAIGHT	10/BOX
JFFS1520S	15.5F X 20CM JET FREE FLOW® ST, STRAIGHT	10/BOX
JFFS1524S	15.5F X 24CM JET FREE FLOW® ST, STRAIGHT	10/BOX
JFFS1528S	15.5F X 28CM JET FREE FLOW® ST, STRAIGHT	5/BOX
JFFS1512IJ	15.5F X 12CM JET FREE FLOW® ST, IJ	5/BOX
JFFS1515IJ	15.5F X 15CM JET FREE FLOW® ST, IJ	5/BOX
JFFS1520IJ	15.5F X 20CM JET FREE FLOW® ST, IJ	5/BOX

**Basic Set Contents:** (1) Catheter with Stylet (1) Guidewire (1) Introducer Needle (3) Dilator (2) End Cap



## IMPORTANT RISK INFORMATION

**Indications For Use:** The Free Flow® ST Catheter is indicated for use in attaining Short-Term vascular access for hemodialysis, hemofiltration, and apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include the subclavian and femoral veins as required. This catheter is indicated for a duration less than (30) days. For femoral placement, monitor catheter condition closely.

**Contraindications:** This catheter is intended for Short-Term vascular access only and should not be used for any purpose other than indicated in these instructions.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit [www.jetmedical.net/patientinformation](http://www.jetmedical.net/patientinformation)



**Medical Components, Inc.**  
1499 Delp Drive, Harleysville, PA 19438, U.S.A.  
P: 215.256.4201 F: 215.256.1787  
[www.medcomp.net](http://www.medcomp.net)



Rue Fritz-Courvoisier 40  
2300 La Chaux-de-Fonds, Switzerland  
T: +41-32-9679021 F: +41-32-\*9679001  
[www.jetmedical.net](http://www.jetmedical.net)

**ISO13485**  
Certified Quality System

PN2740BSI Rev. A 5/24