



JET VALVED PEELABLE INTRODUCER

ACCESSORIES






Wings



Valve



Valve Tip



VALVED PEELABLE INTRODUCER

ACCESSORIES

VALVED PEELABLE INTRODUCERS			
CODE	DESCRIPTION	QTY	GUIDEWIRE ACCEPTANCE
JM15	15F X 14CM JET VALVED PEELABLE INTRODUCER	5/BOX	0.038
JM16	16F X 14CM JET VALVED PEELABLE INTRODUCER	5/BOX	0.038
JM17	17F X 14CM JET VALVED PEELABLE INTRODUCER	5/BOX	0.038

IMPORTANT RISK INFORMATION

Contraindications: The Valved Peelable Introducer is not designed for use in the arterial system or as a hemostatic device.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.jetmedical.net/patientinformation



Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201 F: 215.256.1787
www.medcomp.net



Rue Fritz-Courvoisier 40
2300 La Chaux-de-Fonds, Switzerland
T: +41-32-9679021 F: +41-32-*9679001
www.jetmedical.net

ISO13485
Certified Quality System

PN2752BSI Rev. A 5/24