


JET PERITONEAL DIALYSIS CATHETERS

Accessories



PD Accessories		
JBCAGL	BARBED LUER LOCK ADAPTOR	20/BOX
JTBCGL	TITANIUM BARBED CONNECTOR	50/BOX
JCSGL	PD CATHETER STYLET	5/BOX
JTSGL	PD CATHETER TUNNELER	10/BOX
JFSXGL	X-SERIES CURVED TUNNELER	5/BOX
JTRXGL	X-SERIES TUNNELING ROD	5/BOX

IMPORTANT RISK INFORMATION

I-SERIES & V-SERIES

Indications for Use: The Peritoneal Catheter is indicated for chronic peritoneal dialysis.

Contraindications: This catheter is contraindicated in patients who are not suitable candidates for peritoneal dialysis therapy.

X-SERIES

Indications for Use: The Peritoneal Catheter is indicated for chronic peritoneal dialysis. The extended length of the X-Series Catheter makes it especially applicable for peritoneal dialysis patients when it is necessary to locate the skin exit site remote from the usual lower abdominal region. The catheter may be particularly indicated in patients with obesity, floppy abdominal skin folds, urinary or fecal incontinence, chronic yeast intertrigo, intestinal stomas, or in patients who desire to take deep tub baths.

Contraindications: This catheter is contraindicated in patients who are not suitable candidates for peritoneal dialysis therapy. An upper chest (presternal) exit-site location of the X-Series Catheter may be contraindicated in patients who have had breast augmentation or breast reconstruction or have a tracheostomy. However, this group of patients may be suitable candidates for an upper abdominal exit-site location of the X-Series Catheter.

BARBED LUER LOCK ADAPTOR

Indications for Use: The Medcomp® Barbed Luer Lock Adaptor is indicated for long-term use as a connection point for peritoneal dialysis catheters and male luers.

Contraindications: Use only with catheter provided by Jet Medical S.A. Known or suspected allergies to the device.

PD CATHETER STYLET

Indications for Use: The Jet Peritoneal Catheter Stylet is indicated for transient use to keep a peritoneal dialysis catheter straight and maneuverable during insertion into the peritoneum.

Contraindications: Use only with catheter provided by Jet Medical S.A.

PD CATHETER TUNNELER

Indications for Use: The Jet Peritoneal Catheter Tunneler is indicated for transient use to dissect tissue to create the subcutaneous tunnel opening.


X-SERIES CURVED TUNNELER

Indications for Use: The Jet X-Series Peritoneal Catheter Curved Tunneler is indicated for transient use to dissect tissue to create the subcutaneous tunnel opening.

X-SERIES TUNNELING ROD

Indications for Use: The Jet X-Series Peritoneal Catheter Tunneling Rod is indicated for transient use to dissect tissue to create the subcutaneous tunnel opening.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Jet
Medical^{SA}

Rue Fritz-Courvoisier 40
2300 La Chaux-de-Fonds
Switzerland

T: +41-32-9679021 F: +41-32-*9679001
www.jetmedical.net






Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201 F: 215.256.1787
www.medcomp.net

ISO13485
Certified Quality System

PN2753GL REV. A 4/24

Jet
Medical^{SA}




www.jetmedical.net

I-Series						
Material Silicone		Tip Design Straight & Coiled	French Size 15F	Cuffs Single & Dual		
SIL						
Code	Description	Box Qty	Distal Tip to Cuff (cm)	Cuff to Cuff (cm)	Cuff to Coil (cm)	
CATHETER ONLY STRAIGHT						
JPP41SGL	15F X 41CM JET I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	15	-	-	
JPP46SGL	15F X 46CM JET I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	20	-	-	
JSC31SGL	15F X 31CM JET I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	10.25	-	-	
JSC37SGL	15F X 37CM JET I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	16.25	-	-	
JSC42SGL	15F X 42CM JET I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	21	-	-	
JSC46SGL	15F X 46CM JET I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	25	-	-	
JDC31SGL	15F X 31CM JET I-SERIES CATHETER, DUAL CUFF	5	6.5	3	-	
JDC32SGL	15F X 32.25CM JET I-SERIES CATHETER, DUAL CUFF	5	7	3.5	-	
JDC37SGL	15F X 37CM JET I-SERIES CATHETER, DUAL CUFF	5	10	5	-	
JDC42SGL	15F X 42CM JET I-SERIES CATHETER, DUAL CUFF	5	15	5	-	
JDC47SGL	15F X 47CM JET I-SERIES CATHETER, DUAL CUFF	5	20	5	-	
CATHETER ONLY COILED						
JPP40CGL	15F X 40.25CM JET I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	16.36	-	13	
JPP57CGL	15F X 57.5CM JET I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	13	-	7	
JPP60CGL	15F X 60CM JET I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	15	-	10	
JDC57CGL	15F X 57.5CM JET I-SERIES CATHETER, DUAL CUFF	5	12	6	8	
JDC63CGL	15F X 63CM JET I-SERIES CATHETER, DUAL CUFF	5	19	6	14	
BASIC SET STRAIGHT						
JET-146GL	15F X 46CM JET I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	20	-	-	
JET-146SGL	15F X 46CM JET I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	25	-	-	
JET-242GL	15F X 42CM JET I-SERIES CATHETER, DUAL CUFF	5	15	5	-	
BASIC SET COILED						
JET-160GL	15F X 60CM JET I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	34	-	15	
JET-257GL	15F X 57.5CM JET I-SERIES CATHETER, DUAL CUFF	5	31	6	12	
JET-263GL	15F X 63CM JET I-SERIES CATHETER, DUAL CUFF	5	38.25	6	19	

Catheter Only Contents: (1) Catheter (1) Female Luer (1) Clamp (1) End Cap



Basic Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Guidewire Advancer (1) Syringe (1) Scalpel (1) Tunneler (1) Peelable Introducer (1) Female Luer (1) Clamp (6) Gauze (1) End Cap



V-Series						
Material Silicone		Tip Design Straight & Coiled	French Size 15F	Cuffs Dual		
SIL						
Code	Description	Box Qty	Distal Tip to Cuff (cm)	Cuff to Cuff (cm)	Cuff to Coil (cm)	
CATHETER ONLY STRAIGHT						
JS200101GL	15F X 43CM JET V-SERIES CATHETER, LEFT	5	15	5	-	
JS200102GL	15F X 43CM JET V-SERIES CATHETER, RIGHT	5	15	5	-	
CATHETER ONLY COILED						
JS200201GL	15F X 63CM JET V-SERIES CATHETER, LEFT	5	17.82	5	13	
JS200202GL	15F X 63CM JET V-SERIES CATHETER, RIGHT	5	17.82	5	13	
JS200501GL	15F X 57CM JET V-SERIES CATHETER, LEFT	5	14.10	5	9	
BASIC SET STRAIGHT						
JS200301GL	15F X 43CM JET V-SERIES CATHETER, LEFT	5	15	5	-	
JS200302GL	15F X 43CM JET V-SERIES CATHETER, RIGHT	5	15	5	-	
BASIC SET COILED						
JS200401GL	15F X 63CM JET V-SERIES CATHETER, LEFT	5	17.82	5	13	
JS200402GL	15F X 63CM JET V-SERIES CATHETER, RIGHT	5	17.82	5	13	
JS200502GL	15F X 57CM JET V-SERIES CATHETER, LEFT	5	14.10	5	9	

Catheter Only Contents: (1) Catheter (1) Female Luer (1) Clamp (1) End Cap

Basic Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Guidewire Advancer (1) Syringe (1) Scalpel (1) Tunneler (1) Peelable Introducer (1) Female Luer (1) Clamp (6) Gauze (1) End Cap

X-Series					
Material Silicone		Tip Design Coiled	French Size 15F		
SIL					
Code	Description	Box Qty	Cuff to Coil (cm)–Lower Catheter	Proximal Cuff to End of Catheter (cm)–Upper Catheter	
CATHETER ONLY					
JS20XCGL	15F X 64CM JET X-SERIES CATHETER	5	12.75	28.5	
BASIC SET					
JS20XSGL	15F X 64CM JET X-SERIES CATHETER	5	12.75	28.5	

Catheter Only Contents: (1) Catheter (1) Titanium Barbed Connector (1) Female Luer (1) Clamp (1) End Cap

Basic Set Contents: (1) Catheter (1) Titanium Barbed Connector (1) Introducer Needle (1) Guidewire (1) Guidewire Advancer (1) Syringe (1) Scalpel (1) Tunneler (1) Female Luer (1) Clamp (6) Gauze (1) End Cap (1) Tape Measure