


# JET PERITONEAL DIALYSIS CATHETER DIMENSION GRAPH

MEASUREMENTS ARE REFERENCE DIMENSIONS


## I-Series Straight (Tenckhoff), Single Preperitoneal Cuff



### I-SERIES STRAIGHT (TENCKHOFF), SINGLE PREPERITONEAL CUFF

| Jet Kit           | A<br>Length (cm) | B<br>Tip to Cuff (cm) | C<br>Cuff (cm) | D<br>Cuff to Proximal Tip (cm) | E<br>Sideholes Length (cm) |
|-------------------|------------------|-----------------------|----------------|--------------------------------|----------------------------|
| JPP41S<br>JET-141 | 41               | 15                    | 0.8            | 25                             | 7                          |
| JET-146           | 46               | 20                    | 0.8            | 25                             | 15                         |


## I-Series Straight (Tenckhoff), Single Subcutaneous Cuff



### I-SERIES STRAIGHT (TENCKHOFF), SINGLE SUBCUTANEOUS CUFF

| Jet Kit            | A<br>Length (cm) | B<br>Tip to Cuff (cm) | C<br>Cuff (cm) | D<br>Cuff to Proximal Tip (cm) | E<br>Sideholes Length (cm) |
|--------------------|------------------|-----------------------|----------------|--------------------------------|----------------------------|
| JSC31S             | 31               | 10.25                 | 0.8            | 20                             | 3                          |
| JSC37S             | 37               | 16.25                 | 0.8            | 20                             | 5                          |
| JSC42S             | 42               | 21                    | 0.8            | 20                             | 7                          |
| JSC46S<br>JET-146S | 46               | 25                    | 0.8            | 20                             | 15                         |


## I-Series Straight (Tenckhoff), Dual Cuff



### I-SERIES STRAIGHT (TENCKHOFF), DUAL CUFF

| Jet Kit           | A<br>Length (cm) | B<br>Tip to First Cuff (cm) | C<br>Cuff (cm) | D<br>Cuff to Cuff (cm) | E<br>Cuff to Proximal Tip (cm) | F<br>Sideholes Length (cm) |
|-------------------|------------------|-----------------------------|----------------|------------------------|--------------------------------|----------------------------|
| JDC31S            | 31               | 6.5                         | 0.8            | 3                      | 20                             | 3                          |
| JDC32S            | 32.5             | 7                           | 0.8            | 3.5                    | 20                             | 5                          |
| JDC37S            | 37               | 10                          | 0.8            | 5                      | 20                             | 5                          |
| JDC42S<br>JET-242 | 42               | 15                          | 0.8            | 5                      | 20                             | 7                          |
| JDC47S            | 47               | 20                          | 0.8            | 5                      | 20                             | 15                         |

### I-Series Coiled (Tenckhoff), Single Preperitoneal Cuff



I-SERIES COILED (TENCKHOFF), SINGLE PREPERITONEAL CUFF

| Jet Kit           | A                    | B                       | C                         | D                 | E         | F                         | G               |
|-------------------|----------------------|-------------------------|---------------------------|-------------------|-----------|---------------------------|-----------------|
|                   | Uncoiled Length (cm) | Coiled Tip to Cuff (cm) | Uncoiled Tip to Cuff (cm) | Cuff to Coil (cm) | Cuff (cm) | Cuff to Proximal Tip (cm) | Coiled Tip (cm) |
| JPP40C            | 40                   | 16.36                   | 26                        | 13                | 0.8       | 12                        | 5               |
| JPP57C            | 58                   | 13                      | 31                        | 7                 | 0.8       | 26                        | 5               |
| JPP60C<br>JET-160 | 60                   | 15                      | 34                        | 10                | 0.8       | 25                        | 5               |


### I-Series Coiled (Tenckhoff), Dual Cuff



I-SERIES COILED (TENCKHOFF), DUAL CUFF

| Jet Kit           | A                    | B                             | C                               | D                 | E         | F                 | G                         | H               |
|-------------------|----------------------|-------------------------------|---------------------------------|-------------------|-----------|-------------------|---------------------------|-----------------|
|                   | Uncoiled Length (cm) | Coiled Tip to First Cuff (cm) | Uncoiled Tip to First Cuff (cm) | Cuff to Coil (cm) | Cuff (cm) | Cuff to Cuff (cm) | Cuff to Proximal Tip (cm) | Coiled Tip (cm) |
| JDC57C<br>JET-257 | 57.5                 | 12                            | 31                              | 8                 | 0.8       | 6                 | 18                        | 5               |
| JDC63C<br>JET-263 | 63                   | 19                            | 38.25                           | 14                | 0.8       | 6                 | 17                        | 5               |


### V-Series Straight (Swan Neck Tenckhoff), Dual Cuff



V-SERIES STRAIGHT (SWAN NECK TENCKHOFF), DUAL CUFF

| Jet Kit  | A           | B                      | C         | D                 | E                         | F                     |
|----------|-------------|------------------------|-----------|-------------------|---------------------------|-----------------------|
|          | Length (cm) | Tip to First Cuff (cm) | Cuff (cm) | Cuff to Cuff (cm) | Cuff to Proximal Tip (cm) | Sideholes Length (cm) |
| JS200101 | 43          | 15                     | 0.8       | 5                 | 22                        | 11                    |
| JS200102 |             |                        |           |                   |                           |                       |
| JS200301 |             |                        |           |                   |                           |                       |
| JS200302 |             |                        |           |                   |                           |                       |

### V-Series Coiled (Swan Neck Curl) Dual Cuff



V-SERIES COILED (SWAN NECK CURL) DUAL CUFF

| Jet Kit  | A                    | B                             | C                               | D                 | E         | F                 | G                         | H               |
|----------|----------------------|-------------------------------|---------------------------------|-------------------|-----------|-------------------|---------------------------|-----------------|
|          | Uncoiled Length (cm) | Coiled Tip to First Cuff (cm) | Uncoiled Tip to First Cuff (cm) | Cuff to Coil (cm) | Cuff (cm) | Cuff to Cuff (cm) | Cuff to Proximal Tip (cm) | Coiled Tip (cm) |
| JS200201 | 63                   | 17.82                         | 34                              | 13                | 0.8       | 5                 | 21                        | 5               |
| JS200402 |                      |                               |                                 |                   |           |                   |                           |                 |
| JC200401 |                      |                               |                                 |                   |           |                   |                           |                 |
| JC200402 |                      |                               |                                 |                   |           |                   |                           |                 |